


# Product Information

© 2021 Kuhn Copper Solutions LLC

**The Kuhn all Copper Mask (KCM) and Kuhn all Copper Insert (CMI)** is made of pure copper mesh that sanitizes and filters the air you breathe. The masks and CMI are unique in that they deposit miniscule amounts of copper atoms in the nasal area. **Copper inactivates most corona viruses, flu, and other respiratory viruses and kills bacteria on contact.** Currently, no clinical trails available. Corona virus, flu, respiratory, other viral infections or allergy symptoms may develop, regardless of precautions taken. Read all enclosed instructions before use. For adult use only.


**The Full Face Mask** is ideal for colds, flu, respiratory infections, allergies, and red tide ( Fig 1, 2 )

Fig 1



Full Face Mask

Fig 2



Full Face Mask

Note: The Kuhn all Copper Mask are not intended as surgical, infection control, dental or medical procedure masks and should not be used under any conditions where bodily fluids or blood may be sprayed directly on the mask, as splash may penetrate the mask.

See Page 2 below for Copper Mesh Insert

**Kuhn all Copper Insert (CMI):** Ideal for dentists, first responders, and health care workers. The CMI is a pure copper mesh insert (2.75" x 2.75") designed specifically to be worn beneath a mask, adding a sanitizing component to the conventional mask, converting it to a cidal (virus inactivating) mask and adding an extra layer of protection. Ideal for men with beards who find the full face mask "pinches".

Fig 3



Kuhn all Copper Insert (CMI)

Top arrow: Kuhn all Copper Insert with adhesive mole skin on upper edge.

Bottom arrow: Conventional Mask